

## Aerosol Therapy

# Frog Nebulizer System

Designed to reduce anxiety during treatments, the child-friendly Roscoe Pediatric Frog Nebulizer provides effective treatments with a state-of-the-art piston compressor. Delivers 43 psi at maximum pressure and 16 psi operating pressure (liter flow range 5 - 8 LPM) and weighs only 3.4 pounds. Equipped with a built-in carrying handle and a convenient nebulizer parking post. Comes complete with compressor, nebulizer, child mask, tubing, angled mouthpiece, five extra filters and a five year warranty.



## Product Specifications

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a built-in carrying handle and a convenient nebulizer parking post
- Comes complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
- Five year warranty
- HCPCS Code: E0570\*
- HCPCS Code: A7003\* for nebulizer kits
- HCPCS Code: A7013\* for foam filters

## Technical Specifications

- Power: 120 VAC/60Hz
- Power Consumption: 1.5 amps
- Sound Level: less than 60dBA at 1 meter
- Maximum Compressor Pressure: 29 to 43 psi
- Typical Operating Pressure: 10 to 16 psi
- Typical Operating Flow Rate: 5 ~ 8 LPM
- Medication Capacity: 5 ml
- Particle Size: Less than 5 microns
- Operating Temperature: 50° - 104°F
- Dimensions: 7 5/8" (L) x 6 1/4" (W) x 5 1/3" (H)
- Weight: 3.4 lb.



Carrying Bag and Replacement Filter

Item #	Description	Qty
50000	Roscoe Pediatric Frog Nebulizer System	each
50002	Nylon Carry Bag for Roscoe Frog Compressor System	each
CIF-1128	Foam Filter, 5/8" diameter x 1/2" thick	50/pk

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



**To Order:**  
 (800) 3-ROSCOE  
 fax: (440) 572-4261  
 roscoemedical.com