

Name: \_\_\_\_\_

Date: \_\_\_\_\_

## KANGAROO JOEY FEEDING PUMP INSTRUCTIONS

1. Wash your hands!
2. Hold bag upright and pour in formula. Close cap. (If using the backpack all air must be removed from bag).
3. **It is important that you do not hang the bag more than a couple of inches above the pump.**
4. Remove the plastic protector from the end of the feeding set.
5. Press On/Off key to turn pump on.
6. The screen will now display "LOAD A SET" Open pump door and loop tubing around the pump wheel stretching lightly. The pump will now display "SET LOADED" on the screen.
7. Press and hold prime key. Once all the tubing is filled with formula, release key and press done.
8. Your pump should be preprogrammed and ready to go **however it is very important to confirm all your settings: RATE: \_\_\_\_\_ ML VTBD: \_\_\_\_\_ FLUSH \_\_\_\_\_ HOURS \_\_\_\_\_** (Flush/hours may not apply to you)
9. Press Run. Pump will begin to run until dose is complete.
10. If you need to stop the pump at any time, press the Run hold the run key for 3 seconds.
11. The bags should be rinsed with water between each feeding. Change feeding bags daily.
12. Remember to keep pump plugged in at all times unless feeding away from home. Feeding pump can run for 18 hours on battery power. It will take 8 hours to fully charge the battery.
13. See attachment for instructions on cleaning of pump.

**If you have any questions, please call Sherwood Clinical at  
800-847-3987**