

Name: \_\_\_\_\_

Date: \_\_\_\_\_

## **ZEVEK FEEDING PUMP INSTRUCTIONS**

1. Wash hands!
2. Hold bag upright and pour in formula. Close cap. If using the backpack all air must be removed from bag.
3. Open pump door and loop tubing around the pump wheel stretching lightly.
4. Press and hold On/Off key for 1.5 seconds to turn pump on.
5. You will need to confirm all your settings: RATE: \_\_\_\_\_ ML DOSE: \_\_\_\_\_
6. Press and hold prime key. Once all the tubing is filled with formula, release key.
7. Press Run/Pause key. Pump will begin to run until dose is complete.
8. If you need to stop the pump at any time, press the Run/Pause key.
9. If you need to change the rate, press Rate/Dose Key. Then press + or – until you reached the desired rate. Then press Rate/Dose key again to display dose. Dose is the total amount of formula that will be given at the feedings. If you need to change the dose, press + or – until you reached the desired dose.

**BEFORE CHANGING YOUR RATE OR DOSE ALWAYS MAKE SURE YOU ARE ON THE CORRECT SCREEN.**

10. The bags should be rinsed with water between each feeding. Change feeding bags daily.
11. Remember to keep pump plugged in at all times unless feeding away from home. Feeding pump can run for 24 hours on battery power. It will take 6 hours to fully charge the battery.
12. See attachment for instructions on cleaning of pump.

**If you have any questions, please call Sherwood Clinical at  
800-847-3987**

