

Aerosol Therapy

Roscoe Bunny Pediatric Nebulizer System

Designed to reduce anxiety during treatments, the child-friendly Roscoe Bunny Nebulizer provides effective treatments with a state-of-the-art piston compressor. Delivers 43 psi at maximum pressure and 16 psi operating pressure (liter flow range 5 - 8 LPM) and weighs only 3.4 pounds. Equipped with a built-in carrying handle and a convenient nebulizer parking post. Comes complete with compressor, nebulizer, child mask, tubing, angled mouthpiece, five extra filters and a five year warranty. The Roscoe Bunny Nebulizer Compressor is also available with the TruNeb Reusable Neb Kit option.



Attractive Retail Packaging!



NEB-BUNNY



Carrying Bag and Replacement Filter

Product Specifications

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a convenient nebulizer parking post
- Comes complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
- TruNeb Reusable Neb Kit available - purchase with the compressor or alone
- BUNNY-TRU includes TruNeb Reusable Neb Kit
- Five year warranty
- HCPCS Code: E0570*
- HCPCS Code: A7003* for disposable nebulizer kits
- HCPCS Code: A7005* for reusable nebulizer kits
- HCPCS Code: A7013* for foam filters

Technical Specifications

- Power: 120 VAC/60Hz
- Power Consumption: 1.5 amps
- Sound Level: less than 60dBA at 1 meter
- Maximum Compressor Pressure: 25 to 42 psi
- Typical Operating Pressure: 8 to 13 psi
- Typical Operating Flow Rate: 5 ~ 8 LPM
- Medication Capacity: 5 ml
- Particle Size: Less than 5 microns
- Operating Temperature: 50o - 104oF
- Dimensions: 8.78" (L) x 5.4" (W) x 9" (H)
- Weight: 3.4 lb.

Item #	Description	Qty
NEB-BUNNY	Roscoe Bunny Pediatric Nebulizer System with Disposable Neb Kit	each
BUNNY-TRU	Roscoe Bunny Pediatric Nebulizer System with one Disposable Neb Kit and one TruNeb Reusable Neb Kit	each
BAG-BUNNY	Nylon Carry Bag for Roscoe Bunny Pediatric Compressor System	each
CIF-1128PK6	Foam Filter, 5/8" diameter x 1/2" thick	6/pk

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:

(800) 3-ROSCOE
fax: (440) 572-4261
roscoemedical.com