

Aerosol Therapy

Roscoe Dinosaur Pediatric Nebulizer System

Designed to reduce anxiety during treatments, the child-friendly Roscoe Dinosaur Nebulizer provides effective treatments with a state-of-the-art piston compressor. Delivers 42 psi at maximum pressure and 13 psi operating pressure (liter flow range 5 - 8 LPM) and weighs only 2.85 pounds. Equipped with a built-in carrying handle and a convenient nebulizer parking post. Comes complete with compressor, nebulizer, child mask, tubing, angled mouthpiece, five extra filters and a five year warranty. The Roscoe Dinosaur Nebulizer Compressor is also available with the TruNeb Reusable Neb Kit option.

Product Specifications

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a convenient nebulizer parking post
- Comes complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
- TruNeb Reusable Neb Kit available - purchase with the compressor or alone
- DINO-TRU includes TruNeb Reusable Neb Kit
- Five year warranty
- HCPCS Code: E0570*
- HCPCS Code: A7003* for disposable nebulizer kits
- HCPCS Code: A7005* for reusable nebulizer kits
- HCPCS Code: A7013* for foam filters

Technical Specifications

- Power: 120 VAC/60Hz
- Power Consumption: 1.5 amps
- Sound Level: less than 60dBA at 1 meter
- Maximum Compressor Pressure: 25 to 42 psi
- Typical Operating Pressure: 8 to 13 psi
- Typical Operating Flow Rate: 5 ~ 8 LPM
- Medication Capacity: 5 ml
- Particle Size: Less than 5 microns
- Operating Temperature: 50°F - 104°F
- Dimensions: 4.06" (L) x 11.42" (W) x 5.2" (H)
- Weight: 2.85 lb.



NEB-DINO



Attractive Retail Packaging!

**** Preorder item for December Delivery ****

Item #	Description	Qty
NEB-DINO	Roscoe Dinosaur Pediatric Nebulizer System with Disposable Neb Kit	each
DINO-TRU	Roscoe Dinosaur Pediatric Nebulizer System with one Disposable Neb Kit and one TruNeb Reusable Neb Kit	each
BAG-DINO	Nylon Carry Bag for Roscoe Dinosaur Pediatric Compressor System	each
CIF-1128PK6	Foam Filter, 5/8" diameter x 1/2" thick	6/pk

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com