

# Aerosol Therapy

## Dog Nebulizer System

Designed to reduce anxiety during treatments, the child-friendly Roscoe Dog Nebulizer provides effective treatments with a state-of-the-art piston compressor. Delivers 43 psi at maximum pressure and 16 psi operating pressure (liter flow range 5 - 8 LPM) and weighs only 3.4 pounds. Equipped with a built-in carrying handle and a convenient nebulizer parking post. Comes complete with compressor, nebulizer, child mask, tubing, mouthpiece, five extra filters and a five year warranty.



### Product Specifications

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a built-in carrying handle and a convenient nebulizer parking post
- Comes complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
- Five year warranty
- HCPCS Code: E0570\*
- HCPCS Code: A7003\* for nebulizer kits
- HCPCS Code: A7013\* for foam filters

### Technical Specifications

- Power: 120 VAC/60Hz
- Power Consumption: 1.5 amps
- Sound Level: less than 60dBA at 1 meter
- Maximum Compressor Pressure: 29 to 43 psi
- Typical Operating Pressure: 10 to 16 psi
- Typical Operating Flow Rate: 5 ~ 8 LPM
- Medication Capacity: 5 ml
- Particle Size: Less than 5 microns
- Operating Temperature: 50° - 104°F
- Dimensions: 8 1/8" (L) x 5 1/2" (W) x 9" (H)
- Weight: 3.4 lb.



Carrying Bag and  
Replacement Filter

Item #	Description	Qty
NEB-DOG	Roscoe Pediatric Dog Nebulizer Compressor System	each
BAG-NBDG	Nylon Carry Bag for Roscoe Dog Nebulizer System	each
CIF-1128	Foam Filter, 5/8" diameter x 1/2" thick	50/pk

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:  
(800) 3-ROSCOE  
fax: (440) 572-4261  
roscoemedical.com